

550 Disposable Pipette Tips

Apricot designs and manufactures 550 µl disposable pipette tips made of premium research grade polypropylene. These tips are compatible with the Apricot multichannel pipettors and offer precise pipetting for mid to high volumes.

- Highly hydrophobic -- minimum fluid retention
- Thin wall at tip end -- increases accuracy
- Small size at tip end -- ideal for 384 and 96-well labware
- Sterilization -- autoclavable to 120C
- High quality -- stringent manufacturing standards
- Uniquely packaged -- racks individually packaged for cleanliness and ease of handling

The Apricot 550 tips are available in two different styles of tip packaging. The R style packaging can be used with Apricot EZ-Load heads and the B style package is suitable for robotic handling for high-throughput applications. Both tip packages feature a retention collar that keeps the tips in alignment while they are loaded on the pipettor head.



“R” Style Tip Packaging

Compatible Heads:

096-A-01-EZL-550 (EZ-LOAD)

096-A-01

550-A-01N

125-A-05-55-41DN



“B” Style Tip Packaging

Compatible Heads:

096-A-01

550-A-01N

125-A-05-55-41DN

Ordering Information

Part Number	Description	Max Volume	Tips per Rack	Racks per Case
550-96R-NS	Non-Sterile	550 ul	96	20
550-96R-S	Sterile	550 ul	96	20
550-96R-FS	Filter Tip, Sterile	400 ul	96	20
550-96B-NS	Non-Sterile	550 ul	96	20
550-96B-S	Sterile	550 ul	96	20
550-96B-FS	Filter Tip, Sterile	400 ul	96	20