

Apricot Total Pipetting Solution (TPS) Selection Guide

	TPS-24	TPS-96H	TPS-96L	TPS-384
No. of Labware Stations	6	12,15	12,15	12,15
No. of Pipetting Channels	24	96	96	384
Max Capacity per Channel (ul)	125	550	125	125
96-well	√	√	√	√
384-well	√	√	√	√
1536-well	√			√
Serial Dilution col to col	√	√	√	√
Serial Dilution row to row	√	√	√	√

Apricot Personal Pipettors™ (PP) Multistation Selection Guide

	PP-96H-M	PP-96L-M	PP-96H-H	PP-96S-100	PP-384-M
No. of Labware Stations	6	6	6	6	6
No. of Pipetting Channels	96	96	96	96	384
Max Capacity per Channel (ul)	550	125	550	100	125
96-well	√	√	√	√	√
384-well	√	√	√	√	√
1536-well					√
Serial Dilution col to col	√	√	√		√
Serial Dilution row to row	√	√	√		√
SPE			√		

Apricot i-Pipette Single Station Selection Guide

	i-Pipette 96-500	i-Pipette 96-125	i-Pipette 384	i-Pipette Pro 96-500	i-Pipette Pro 96-125	i-Pipette Pro 384
No. of Labware Stations	1 or 2	1 or 2	1 or 2	1 or 2	1 or 2	1 or 2
No. of Pipetting Channels	96	96	384	96	96	384
Max Capacity per Channel (ul)	500	125	125	500	125	125
96-well	√	√	√	√	√	√
384-well	√	√	√	√	√	√
1536-well						
Serial Dilution col to col	√	√	√	√	√	√
Serial Dilution row to row	√	√	√	√	√	√
Motorized Head Clamp and Motorized Elevator				√	√	√

Designs that make sense

Apricot Designs, Inc. 681 Arrow Grand Circle, Covina, CA 91722, USA Tel: (626) 966-3299

