



The most innovative programmable pressure processor available in a stand-alone unit.

Delivering precise, uniform flow control across multiple ports for processing 96-well extraction plates or individual columns.

## Flexibility · Control · Reproducibility

SPE/SLE Positive Pressure Processor

### Why choose the EVEN96?

Proven Engineering Expertise from Apricot Designs. Helping you get the best results for over 25 years.

### Innovative Design

Manufactured with flexibility, control and reproducibility in mind.

### Intuitive and Accurate Control

Interface Control at your fingertips via Microsoft Surface® Control Tablet.  
No external computer needed!

### Flexible Processing of Plates/Columns

Readily process 96-well format SPE/SLE plates as well as individual SPE/SLE columns on the same platform.

### Precise Height Control & Superior Sealing

Engineered for incremental control and accurate positioning.

### Precision Flow Control

Optimized programmable flow controller. Provides a flow range from microliters per minute ( $\mu\text{L}/\text{min}$ ) for flow-sensitive methods to liters per minute L/min for sorbet drying.

### High Performance Internal Heater

Internal heater for fast sorbent drying in complex extraction protocols.

### Advanced Safety Features

Manifold pinch sensor for enhanced operator safety. Fits within suitably sized laboratory hoods for solvent vapor control.