



The most innovative programmable pressure processor available in a stand-alone unit.

Delivering precise, uniform flow control across multiple ports for processing 96-well extraction plates or individual columns.

Flexibility · Control · Reproducibility

SPE/SLE Positive Pressure Processor

Why choose the EVEN96?

Proven Engineering Expertise from Apricot Designs. Helping you get the best results for over 25 years.

Innovative Design

Manufactured with flexibility, control and reproducibility in mind.

Intuitive and Accurate Control

Interface Control at your fingertips via Microsoft Surface® Control Tablet.
No external computer needed!

Flexible Processing of Plates/Columns

Readily process 96-well format SPE/SLE plates as well as individual SPE/SLE columns on the same platform.

Precise Height Control & Superior Sealing

Engineered for incremental control and accurate positioning.

Precision Flow Control

Optimized programmable flow controller. Provides a flow range from microliters per minute ($\mu\text{L}/\text{min}$) for flow-sensitive methods to liters per minute L/min for sorbet drying.

High Performance Internal Heater

Internal heater for fast sorbent drying in complex extraction protocols.

Advanced Safety Features

Manifold pinch sensor for enhanced operator safety. Fits within suitably sized laboratory hoods for solvent vapor control.