



Ask about our  
30th Anniversary ~  
Limited Edition  
3-station w/iPad mini®

# iPIPETTEPRO

Proven high-performance automated liquid handling

High-performance liquid handling  
EZ-Load Tip Technology  
Automated X & Z positioning  
Small footprint

## FEATURES

3 Stations  
30th Anniversary Badge  
Quick Inter-changeable Head



## FEATURES

- Proven high-performance automated liquid handling – 8,12,16, 24, 96, 384 channel pipetting
- Control via easy to use interface – Apple iPod® or Apple iPad mini® controller
- EZ-Load® Tip Technology – Load tips quickly and easily, without tools!
- Automated plate loading via motorized elevator
- 2 or 3-position shuttle /programmable serial dilution

## MODELS

- 384 0.5–125µL
- 96H 1–500µL
- 96L 0.5–125µL
- 96E 5–1000µL

*Ask us about custom configurations and accessories*

## APPLICATIONS

- Plate replication
- Plate reformatting
- Reagent addition
- Compound addition
- Serial dilution by column

## SPECIFICATIONS

- Weight: 45–52lbs
- Width: 9”–15” Depth: 14” Height: 18.5”– 20”
- Power 100–240VAC, 50/60Hz
- Operating temperature 50–100°F

**[www.apricotdesigns.com](http://www.apricotdesigns.com)**

Apricot Designs, Inc. 677 Arrow Grand Circle, Covina, CA 91722 USA

Telephone: 626-966-3299 • Fax: 626-966-3200

**For product information, quotes, or any other questions, feel free to e-mail us at [info@apricotdesigns.com](mailto:info@apricotdesigns.com)**

iPod® or Apple iPad mini® are registered trademarks of Apple Inc. The Apricot Designs logo is a trademark of Apricot Designs, Inc. © 2019 Apricot Designs, Inc. Specifications subject to change. All information is deemed accurate at the time of publication. Apricot Designs is not responsible for any errors that may occur.