

# PERSONAL PIPETTOR +Walkaway a work cell resource

Easy to use high-performance pipettor



**Apricot Designs' Personal Pipettor +Walkaway** combines the functionality of several pipettors into a single work cell resource. Complete multi-step assays of non-similar tip formatting with a single click.

The Personal Pipettor +Walkaway is a continuous automated work cell with an intuitive operator interface and automated tip changes to help you get the best results - accurate, reliable, and repeatable.

The Personal Pipettor +Walkaway minimizes repetitive tasks, like serial dilutions, eliminates manual pipetting errors, and delivers results in a small footprint that maximizes bench space.

## PIPETTING JUST GOT BETTER

### APPLICATIONS

- Plate replication
- Plate reformatting
- Reagent addition
- Compound addition
- Serial dilution by row or column
- Purification
- PCR Prep
- ELISA
- NGS library preparation
- Cell-based assay
- +Much more! ...

### FEATURES

- 5 linear stations, plus 1 tip wash/reservoir station
- 1, 8, 12, 16, 24, 96, 384 channel pipetting
- Multiple pipetting heads of various formats on deck
- High-precision pipetting heads are independent of pipetting mechanism
- On the fly tip changes with EZ-Load tip technology
- Intuitive iPad® controller is easy to use. Also, available with PC software controller
- Microplate stacker and robotics ready - readily integrates into an automation work cell
- Integrated tip washing system
- Small footprint maximizes bench space and fits in most standard hoods

### SPECIFICATIONS

Dimensions	28"W x 19"D x 24"H 711mmW x 483mmD x 607mmH
Weight Approx.	98lbs (45kg)
Power	100-240 VAC, 50/60 Hz
Operating Temp	50-100°F (10-37°C)

[www.apricotdesigns.com](http://www.apricotdesigns.com)

Apricot Designs, Inc. 677 Arrow Grand Circle, Covina, CA 91722 USA  
Telephone: 626-966-3299 | Fax: 626-966-3200

Apricot Designs (Shanghai) Co., Ltd. 上海美杏高德生物科技有限公司  
上海市浦东新区张江路 665号德宏大厦1008室 邮编:201203

电话:(86)21 50796338  
传真:(86)21 50796228

© 2018 Apricot Designs, Inc. The Apricot Designs logo is a trademark of Apricot Designs, Inc.  
Specifications subject to change. All information is deemed accurate at the time of publication.  
Apricot Designs is not responsible for any errors that may occur.