



An innovative new product bringing automated liquid handling within reach of labs that performs pipetting operations onto microplates.

Innovation · Accuracy · Precision

Why choose the S-Pipette?

Proven automated pipetting engineering from Apricot Designs.
Helping you get the best results for over 25 years.

Innovative design

Manufactured with flexibility, versatility and affordability in mind

Intuitive and Accurate Control Interfaces

Control at your fingertips, or in the palm of your hand, through a variety of interfaces:

- Touch screen ■ Quick access keys
- Apple® ipod or iPhone

No external computer needed!

Interchangeable Pipetting Core

Use the 96 or 384 channel pipetting configuration that best matches the applications you perform

Adjustable Split Level Plate Elevators

Control individual sampling using multiple plates and a wide variety of liquid handling labware
2, 3 or 4 Plate Elevators available.

Additional Flexibility, manual indexing of the pipetting head

Access 384-well plates with a 96-channel version or access 96-well plates with a 384-channel version (with 96-channel head)

Applications

Plate Replication • Plate Reformatting
Reagent Addition • Compound Addition
Serial Dilution • Indexing by Column